





Before testing, scan the QR code to watch the "how to use" video.

## BEFORE STARTING

Wash your hands with soap and water for at least 20 seconds before and after test. If soap and water are not available, use hand sanitizer with at least 60% alcohol.



## 1. PREPARE FOR THE TEST

### 1A. Check the expiration date on the box.

Do not use if the kit if it has been damaged or has expired.

### 1B. Ensure kit is at room temperature for at least 30 minutes prior to use.

Open the box carefully as it will be used in a later step (1D). Do not open individual components until instructed.

**Note:** A timing device (clock, timer, phone etc.) is required, but not provided.

### 1C. Remove the cover of the tube with Extraction buffer.




### 1D. Put the Tube in the tube holder in the box.

**Note:** Being careful not to spill the Tube contents.

## 2. NASAL SWAB SPECIMEN COLLECTION

### 2A. Open Swab protective pouch.


Remove the sterile swab from the pouch



**!** Keep fingers away from the Swab end. Touch the stick end only.

### 2B. Swabbing both nostrils.

Insert the soft end of the Swab into your nostril until you feel resistance (Approx. 2cm up your nose).



Slowly twist the swab, rubbing it along the insides of your nostril, 5 - 10 times against the nasal wall. Gently remove Swab from nostril.

### 2C. Using the same Swab, repeat step 2B, in your other nostril.

Withdraw the swab.


#### Note:

- 1 This may feel uncomfortable. Do not insert the swab any deeper if you feel strong resistance or pain.
- 2 When the nasal mucosa is damaged or bleeding, nasal swab collection is not recommended.
- 3 If you are swabbing others, please wear a face mask. With children, you may not need to insert the swab as far into the nostril.
- 4 For very young children, you may need another person to steady the child's head while swabbing.

### 2D. Insert the Swab into the extraction Tube.

Ensure it is touching the bottom and stir the swab to mix well.

Press the swab head against the tube and rotate the swab for 10 - 15 seconds




### 2E. Hold the Tube firmly with one hand.

Remove the swab while squeezing the swab head against the inside of the Extraction tube. Place the swab in a plastic bag.

## 2F. Close the cap of the extraction tube

Return the Tube to the Kit Box tube holder before proceeding to the next step.



## 3. PERFORM THE TEST


### 3A. Remove the test cassette from the sealed foil pouch and use it within one (1) hour.

**Note:** Best results will be obtained if the test is performed immediately after opening the foil pouch.

Place the test cassette on a flat and level surface.

**!** Do not move the test cassette during test developing.

### 3B. Invert the specimen extraction tube and add 3 drops of extracted specimen to both sample well (S) of the test cassette.




Start the timer. Secure tube cap back on extraction tube and wait 10 minutes.

**!** Do not touch the Test Device during this period.

### 3C. Read the result at 10 minutes.

Keep Test Device flat on table. Do not read the result earlier than 10 minutes or after 20 minutes.



## 4. READING THE RESULTS

Please share your test result with your healthcare provider and carefully follow your local COVID guidelines/requirements.

### POSITIVE

COVID-19



**POSITIVE SARS-CoV-2:\* Two colored lines appear in the COVID-19 window.**

One colored line should be in the control region (C) and another colored line should be in the Test region (T).

FLU A+B



**POSITIVE Influenza A:\* Two colored lines appear in the FLU A+B window.**

One colored line should be in the control region (C) and another colored line should be in the Influenza A region (A).

FLU A+B



**POSITIVE Influenza B:\* Two colored lines appear in the FLU A+B window.**

One colored line should be in the control region (C) and another colored line should be in the Influenza B region (B).


FLU A+B



**POSITIVE Influenza A and Influenza B:\* Three colored lines appear in the FLU A+B window.**

One colored line should be in the control region (C) and two colored line should be in the Influenza A region (A) and Influenza B region (B).

COVID-19 FLU A+B



**POSITIVE SARS-CoV-2 and Influenza A:\* Two colored lines appear in the COVID-19 window; and two colored lines appear in the FLU A+B window.**

COVID-19 window: One colored line should be in the control region (C) and another colored line should be in the Test region (T).

FLU A+B window: One colored line should be in the control region (C) and another colored line should be in the Influenza A region (A).

COVID-19 FLU A+B



**POSITIVE SARS-CoV-2 and Influenza B:\* Two colored lines appear in the COVID-19 window; and two colored lines appear in the FLU A+B window.**

COVID-19 window: One colored line should be in the control region (C) and another colored line should be in the Test region (T).

FLU A+B window: One colored line should be in the control region (C) and another colored line should be in the Influenza B region (B).

COVID-19 FLU A+B



**POSITIVE SARS-CoV-2 and Influenza A+B:\* Two colored lines appear in the COVID-19 window; Three colored lines appear in the FLU A+B window.**

COVID-19 window: One colored line should be in the control region (C) and another colored line should be in the Test region (T).

FLU A+B window: One colored line should be in the control region (C) and two colored line should be in the Influenza A region (A) and Influenza B region (B).

**\*NOTE:** The intensity of the color in the test line region (T/B/A) will vary based on the amount of SARS-CoV-2 and/or Influenza A+B antigen present in the sample. So any shade of color in the test region (T/B/A) should be considered positive.

A positive result means it is very likely you have COVID-19 and/or Influenza A/Influenza B, but the positive samples should be confirmed to reflect this.

POSITIVE SARS-CoV-2: Follow the guidance from your local State or Territory Health Department for guidance on confirmation testing if necessary, and if unwell seek medical assistance

POSITIVE Influenza A and/or Influenza B: Individuals with a positive result or who are unwell are must consult a medical practitioner for follow-up clinical care.

### NEGATIVE

COVID-19 FLU A+B



**NEGATIVE: One colored line appears in the control region (C).**

**No apparent colored line appears in the test line region (T/B/A).**

You are unlikely to have COVID-19 and/or Influenza A/Influenza B. However, it is possible for this test to give a negative result that is incorrect (a false negative) in some people with COVID-19 and/or Influenza A/Influenza B. This means you could possibly still have COVID-19 and/or Influenza A/Influenza B even though the test is negative.

In addition, you can repeat the test with a new test kit. In case of suspicion, repeat the test after 1-2 days, as the coronavirus/Influenza virus cannot be precisely detected in all phases of an infection. If symptoms persist advice to conduct repeat testing and consult a medical practitioner for follow-up clinical care.

Even with a negative test result, distance and hygiene rules must be observed, migration/traveling, attending events and etc. should follow your local COVID/Influenza guidelines/requirements.

### INVALID

COVID-19



FLU A+B



**INVALID: Control line fails to appear.**

Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test or contact with a COVID-19 and/or Influenza test center. If invalid result continues after repeating, advice to contact the sponsor.

## 5. DISPOSE THE TEST KIT

After the test is complete, place all the components in a plastic bag and tightly sealed, then dispose in household waste or rubbish bin.

